

# 3 PLY SURGICAL FACEMASK WITH EAR LOOP

>98% 3-PLY  
NON-STERILE  
BFE FILTRATION

## CST-MASK5098 | CLASS 1 TYPE IIR - Level 3

Certified International EN14683:2019 + AC:2019

Use the masks enclosed to keep your coughs and sneezes covered and to reduce the risk of inhaling airborne contaminants such as infected droplets.

### Instructions:

- Place both straps over the ears
- Mould mask over nose
- Pull mask down to cover chin

Discard mask after use or immediately if it becomes wet or soiled. (Do not re-use)  
Clean your hands immediately after touching used mask.

### Contents:

Adult sized 3ply disposable masks– single use

Certified International Standard EN14683:2019+ AC:2019

Tested to >98% bacterial filtration efficiency.

Masks are one size fits all but can be moulded to fit bridge of nose .



TEST	TEST STANDARD	STANDARD	TEST DATE	VALUE
Bacterial Filtration Efficiency: BFE	EN14683:2019 (AnnexeB)	>98% (Type IIR)	June, 18 2020	99.9%min
Breathability (DELTA P)	EN14683:2019 (AnnexeC) +AC: 2019	<60Pa/cm <sup>2</sup> (Type IIR)	June, 18 2020	48max
Resistance	ISO22609:2004	≥16kPa	June, 18 2020	Compliant
Intracutaneous irritation test	ISO10993-10	no irritation	June, 18 2020	No irritation
Cytotoxicity	ISO10993-5	No cytotoxicity	June, 18 2020	No cytotoxicity
Sensitization	ISO10993-10 et-10A	Topical application: no sensitization	June, 18 2020	Topical application: no sensitization
Sensitization	ISO10993-10 et-10A	Intradermal injection: no sensitization	June, 18 2020	Intradermal injection: no sensitization
Microbial cleanliness CFU/g	NF11737:2018	≤30 (Type IIR)	June, 18 2020	13-18 CFU/g
Using time: BFE+Delta P	EN14683:2014 (AnnexeB&C)	>98%(Type IIR) <49Pa/cm <sup>2</sup>	June, 18 2020	8h:>99.9% 8h:<37,00Pa/cm <sup>2</sup>

## TECHNICAL INFORMATION

Usual Colour	Blue
Outer Layer Component	100% Blue SBPP-23gsm
Filtrate Component	PP White Filtrate Paper (Meltblown) - 22 gsm
Inner Layer Component	100% White SBPP-23gsm
Elastic	No Latex round white elastic
Mask Length	17.5cm
Mask width	Before Pleat - 17cm After Pleat - 9.5 cm
Aluminium Length	10.5cm
Relaxed Ear Loop	16cm
Ear Loop Elastic Width	0.6cm

## PACKING

Inner Packing	50 pieces sealed in an inner box
Outer Packing (shipper)	40 boxes/shipper (2000 masks)
Shipper Size/Weight	52 x 38 x 22cm - 7kg
Part Number	CST- Mask5098



NATIONAL DISTRIBUTION

Phone: +61 425 236 248

Email: [info@cleanstreamtechnologies.com](mailto:info@cleanstreamtechnologies.com)

Visit: [www.cleanstreamtechnologies.com](http://www.cleanstreamtechnologies.com)



CLEAN STREAM TECHNOLOGIES