







Product Description		Features					Measurements			Physical Properties			Packaging		Standards		
Glove Name	Product Codes	Material	Colour	Texture	Powder Free	Hand Specific	Length (mm)	Thickness (mm)			Force at break (N) Before Ageing	Elongation at break (%) Before Ageing	Modulus at 300% elongation (N) Before Ageing	Box Qty	Carton Qty	Chemo Tested	Australian Surgical Glove Standards Compliance
								Cuff	Palm	Finger							
HAMILTON 	SNL2855 - 5.5 SNL2860 - 6.0 SNL2865 - 6.5 SNL2870 - 7.0 SNL2875 - 7.5 SNL2880 - 8.0 SNL2885 - 8.5 SNL2890 - 9.0	Latex	Natural White 	Fully	✓	Hand Specific (Curved Fingers)	290 ± 10	0.15 ± 0.04	0.19 ± 0.04	0.21 ± 0.04	≥ 12.5	≥ 700	≥ 2.0	50 pairs (5.5-9.0)	4	✓	AS/NZS 4179:2014 ASTM D3577 EN 455 part 1, 2, 3 & 4
VICTOR OUTER 	SPI2855 - 5.5 SPI2860 - 6.0 SPI2865 - 6.5 SPI2870 - 7.0 SPI2875 - 7.5 SPI2880 - 8.0 SPI2885 - 8.5 SPI2890 - 9.0	Polyisoprene	Natural White 	Fully	✓	Hand Specific (Curved Fingers)	290 ± 10	0.15 ± 0.04	0.19 ± 0.04	0.21 ± 0.04	≥ 12.5	≥ 700	≥ 2.0	50 pairs (5.5-9.0)	4	✓	AS/NZS 4179:2014 ASTM D3577 EN 455 part 1, 2, 3 & 4
VICTOR UNDER 	SPI2255 - 5.5 SPI2260 - 6.0 SPI2265 - 6.5 SPI2270 - 7.0 SPI2275 - 7.5 SPI2280 - 8.0 SPI2285 - 8.5 SPI2290 - 9.0	Polyisoprene	Blue 	Fully	✓	Hand Specific (Curved Fingers)	290 ± 10	0.15 ± 0.04	0.19 ± 0.04	0.21 ± 0.04	≥ 12.5	≥ 700	≥ 2.0	50 pairs (5.5-9.0)	4	✓	AS/NZS 4179:2014 ASTM D3577 EN 455 part 1, 2, 3 & 4